Instructions for Use

DownStream[®] System DS-2 DownStream[®] Cartridge DSC-2EU

- Read all instructions and warnings prior to use.
- The DownStream[®] System equipment may be used only by trained personnel under the supervision of physicians trained in angiography and percutaneous coronary intervention (PCI).
- Do not re-sterilize or reuse the DownStream[®] Cartridge, SSO₂ delivery catheter, or disposable accessories; these devices are for single use only.
- Only TherOx personnel or authorized representatives may ship, install and service the DownStream System Console "Console".

Device Description

SuperSaturated Oxygen (SSO₂) Therapy is an adjunctive cardiac catheterization laboratory initiated procedure targeted at the left main coronary artery (LMCA) of an acute myocardial infarction (AMI) patient after successful percutaneous intervention (PCI) with stenting has been performed of the left anterior descending coronary artery. The equipment necessary for SSO₂ Therapy is comprised of three components: the DownStream[®] System Console (“console”), the DownStream[®] Cartridge (“cartridge”), and the SSO₂ delivery catheter. The console and cartridge function together to create a highly oxygen-enriched saline solution called SuperSaturated Oxygen (“SSO₂”) solution. A small amount of autologous blood is mixed with the SSO₂ solution producing oxygen-enriched hyperoxemic blood which is then delivered to the targeted major epicardial artery via the SSO₂ delivery catheter. The duration of SSO₂ Therapy is 60 minutes.

DownStream Console

The console is the electromechanical device that controls the cartridge and monitors performance and safety during administration of SSO₂ Therapy. The console has safety features that continuously monitor system parameters such as the blood flow rate and pressure and detect potentially unsafe conditions such as the presence of air-in-line. A display screen guides the health care professional through setup and clinical operation. The console is non-sterile and has no direct contact with the patient or the blood flow path. The console is intended to be mains-operated (AC-powered) and stationary during use but is equipped with battery backup power. The console is designed for operation in a cardiac catheterization laboratory and must be operated and monitored by trained personnel when providing SSO₂ Therapy.

DownStream Cartridge

The cartridge is a single-use disposable device that is loaded into the console by a trained healthcare professional. The cartridge has a three-chambered body that creates SSO₂ solution from inputs of hospital-supplied oxygen and physiologic saline and mixes the SSO₂ solution with arterial blood within the cartridge blood path. The cartridge has a tube set that draws the patient’s arterial blood through the draw line and returns oxygen-enriched hyperoxemic blood through the return line to the SSO₂ delivery catheter. The cartridge draw line connects to an arterial sheath. Sheath placement may be coaxial (single arterial access site) or contralateral (two arterial access sites) at the physician’s discretion. A physician makes two line connections during the initiation of SSO₂ Therapy: the cartridge draw line is attached to the arterial sheath before priming the blood flow path, and the return line is attached to the SSO₂ delivery catheter after the blood flow path is successfully primed. The priming volume of the cartridge is approximately 60 ml.

SSO₂ delivery catheter

The SSO₂ delivery catheter is a 5F (O.D.) over-the-wire catheter that is equipped with a standard luer fitting at the proximal end for attachment to the return line of the cartridge. The SSO₂ delivery catheter has a length of 100 cm and is shaped to facilitate placement in the left main coronary ostium using the PCI arterial access site. The SSO₂ delivery catheter is placed in the ostium of the LMCA using a guidewire by the trained physician and is connected to the cartridge after blood priming.

Only a catheter that is qualified by TherOx may be used as the SSO₂ delivery catheter. The only qualified SSO₂ delivery catheter is the 5F Boston Scientific Impulse[®] angiographic catheter.

The DownStream® System Console, DownStream® Cartridge and SSO₂ delivery catheter Set-Up are shown in **Figures 1 through 5**:

Figure 1. TherOx® DownStream® System Console


Figure 2. Close-up of DownStream® System Console Components


Figure 3. TherOx® DownStream® Cartridge


Figure 4. SSO₂ delivery catheter Set-Up (Coaxial Femoral Access Site Approach)


Figure 5. SSO₂ delivery catheter Set-Up (Radial Delivery Femoral Draw Approach)

Indications for Use

The TherOx DownStream System is indicated for the preparation and delivery of SuperSaturated Oxygen Therapy (SSO₂ Therapy) to targeted ischemic regions perfused by the patient's left anterior descending coronary artery immediately following revascularization by means of percutaneous coronary intervention (PCI) with stenting that has been completed within 6 hours after the onset of anterior acute myocardial infarction (AMI) symptoms caused by a left anterior descending artery infarct lesion.

Contraindications

- Ipsilateral insertion of a second sheath in a single femoral artery for SuperSaturated Oxygen Therapy is strictly contraindicated.
- Presence of an intra-aortic balloon pump.
- Proximal coronary stenosis that restricts flow with the SSO₂ delivery catheter in place.
- Presence of a post-intervention non-stented coronary dissection or perforation.
- Cardiac valvular stenosis or insufficiency, pericardial disease or non-ischemic cardiomyopathy.
- Pregnant or nursing women.
- Cardiogenic shock.
- Patients contraindicated for anticoagulation therapy.
- Subjects with ventricular pseudoaneurysm, VSD, or severe mitral valve regurgitation (with or without papillary muscle rupture).
- Hemoglobin < 10 g/dL.
- Gastrointestinal or genitourinary bleeding within the last two months, or any major surgery (including CABG) within six weeks of procedure.

Summary of Warnings

- Only physicians who have received appropriate training in the use of the DownStream System should perform SuperSaturated Oxygen Therapy (SSO₂ Therapy).
- SuperSaturated Oxygen Therapy must be performed in the cardiac catheterization laboratory (cath lab). The patient must not be moved while the catheter is in place.
- Under no circumstances should the extracorporeal circuit be disconnected and restarted outside the cath lab.
- Do not use any type of catheter for SSO₂ Therapy other than the qualified SSO₂ delivery catheter. Inspect the SSO₂ delivery catheter and packaging for damage and kinking prior to use; do not use if the catheter or packaging is damaged.
- The DownStream System may not be operated concomitantly with magnetic resonance imaging.
- Use aseptic technique throughout procedure.
- Provide adequate anticoagulation prior to and during therapy per clinical practice standards.
- Inspect the cartridge and packaging for damage prior to use; do not use if the cartridge or packaging is damaged.
- Use only 0.9% Normal saline (Isotonic solution).
- Verify that all connections are secure and tight prior to initiating prime.
- Verify that the return line of the cartridge is not connected to the SSO₂ delivery catheter prior to initiating prime.
- The bubble detector is disabled during prime. Ensure that return line is fully blood primed with no presence of air or bubbles in the line prior to making wet-to-wet connection to the SSO₂ delivery catheter.
- Do not disconnect return line from SSO₂ delivery catheter or draw line from sheath when therapy is in progress. Procedure must be stopped before disconnection.
- Do not exceed the SSO₂ Therapy duration of 60 minutes.
- Do not re-use the single-use cartridge and SSO₂ delivery catheter. Re-use of the cartridge and SSO₂ delivery catheter may result in patient infection. In addition, the cartridge is

programmed for single use only and attempted re-use will result in failure to provide treatment.

- Dispose of used cartridges and SSO₂ delivery catheters using standard hospital procedures to eliminate hazardous waste.
- Only qualified personnel trained in oxygen safety and handling should change the oxygen bottle.
- To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
- No modification of this equipment is allowed.

Summary of Precautions

- The index PCI procedure may be done using femoral or radial access; the SSO₂ Therapy blood draw requires a femoral arterial access site.
- If the recommended 7F Merit Medical Custom Sheath Introducer is not used, or less than an 8F introducer sheath (coaxial approach) or 5F introducer sheath (dual site approach) is used for femoral access, SSO₂ Therapy may stop due to inadequate flow rate.
- Do not use the console if the oxygen bottle pressure is less than 800 psig (54.4 bar).
- Do not expose electrical connections, including the transducer port, to fluid contact.
- Do not stretch the return line while loading into the flow probe.
- Ensure that the return line is fully seated in the flow probe slot before latching to prevent the flow probe door from pinching the line.
- Ensure fluoroscopic contrast agents delivered through the SSO₂ delivery catheter have been flushed prior to priming circuit; these viscous solutions may result in a flow stoppage after connection to the cartridge return line is made.
- Ensure that draw and return lines are not kinked or restricted prior to starting prime or at any other time during use.
- Do not unplug the transducer once the cartridge has been prepped.
- After placement, inspect the SSO₂ delivery catheter under fluoroscopic visualization to ensure that it is seated correctly in the LMCA and does not restrict blood flow.
- If the SYSTEM BATTERY POWER LOW message is displayed, plug the console into an electrical outlet.
- Position the console as necessary to provide access to the AC power cable and the Potential Equalization Terminal located on the instrument's back panel. Potential equalization can be managed in accordance with IEC 60601-1.

Potential Adverse Events

Potential adverse events include risks associated with interventional cardiology procedures. Complications from these events may necessitate emergency CABG or intra-aortic balloon pump placement, or re-catheterization.

- Death
- Acute Myocardial Infarction (AMI)
- Stent Thrombosis
- Revascularization (CABG or PCI)
- Congestive Heart Failure
- Hemorrhage
- Abrupt vessel closure/spasm
- Allergic reactions
- Aneurysm
- Anxiety/Dizziness
- Arrhythmias
- Arteriovenous Fistula/Pseudoaneurysm
- Blood loss/damage
- Cardiogenic Shock
- Chest pain/angina

- Coronary Artery Occlusion
- Embolism (including air emboli and thromboemboli)
- Hematoma
- Hemolysis
- Hypertension/Hypotension
- Infection
- Myocardial Rupture
- Nausea/Vomiting
- Neck/back/groin pain
- Pericardial effusion
- Pulmonary Edema
- Renal complications
- Respiratory complications
- Restenosis
- Stroke/TIA
- Tamponade
- Thrombosis
- Vascular damage (dissection, perforation, rupture, or other mechanical injury)

Patient Selection Criteria

- Patient must be ≥ 18 years and < 80 years of age.
- Within 6 hours from symptom onset to successful reperfusion by means of PCI with stenting of left anterior descending coronary artery ST-elevation AMI.
- Successful revascularization by means of PCI with stenting as documented by $< 50\%$ diameter residual stenosis and \geq TIMI flow grade II in the target vessel.
- Systemic arterial pO₂ greater than or equal to 10.7 kPa or 80 mmHg with supplemental oxygen.

Equipment Requirements

Quantity	Equipment
1	DownStream System Console– Model DS-2
1	DownStream Cartridge – Model DSC-2EU
1	SSO ₂ delivery catheter*
1 liter	Physician-prescribed, sterile physiologic 0.9% normal saline
1	"E" cylinder, or equivalent size of in-hospital, medical-grade oxygen with at least 800 psig (54.4 bar) pressure.
1	Guidewire, maximum OD 0.038 inch (0.97mm)
1	Coaxial approach: 7F** Merit Medical Custom Sheath Introducer Kit (K15-00147) or 8F introducer sheath with sidearm- or - Contralateral approach: 5F (minimum) introducer sheath
7-8	3 ml syringes
1-2	5 ml syringes

* 5F Boston Scientific Impulse® angiographic catheter

** 7F Merit Medical Custom Sheath Introducer Kit has a larger diameter sidearm to allow for required high flow blood draw. No other 7F sheath can be used.

Directions For Use

Patient Anticoagulation

WARNING: Provide adequate anticoagulation prior to and during therapy per clinical practice standards.

Patient anticoagulation throughout the SSO₂ Therapy procedure is aligned with current practice in STEMI patients treated with primary PCI as outlined below. Patient anticoagulation may consist of one of the following regimens:

- Heparin and routine use of GP IIb/IIIa inhibitors
- Heparin and cangrelor, with provisional (bail-out) use of GP IIb/IIIa inhibitors allowed
- Bivalirudin and cangrelor, with provisional (bail-out) use of GP IIb/IIIa inhibitors allowed
- Bivalirudin, with provisional (bail-out) use of GP IIb/IIIa inhibitors allowed

Procedural anticoagulation may consist of either heparin with either cangrelor or with routine use of a GP IIb/IIIa inhibitor, or bivalirudin with or without the addition of cangrelor (in which case GPIIb/IIIa inhibition is not recommended unless required for refractory procedural thrombotic complications). If a GP IIb/IIIa inhibitor is used, it should be continued during the SSO₂ infusion and for a total of at least 12 hours post procedure, as per standard of care. Intravenous cangrelor may also be used according to label per physician discretion, administered as a bolus plus infusion, with the infusion continued for 2-4 hours post-PCI.

If unfractionated heparin is selected, local dosing regimens may be used, although it is recommended that an initial bolus of 60 IU/kg is administered, followed by titration to the ACT of 200-250 seconds in patients receiving a GP IIb/IIIa inhibitor.

If bivalirudin is selected as the procedural anticoagulant, it should be administered as a bolus of 0.75 mg/kg IV prior to PCI, followed by an infusion of 1.75 mg/kg/h initiated as soon as possible. Intravenous cangrelor may also be used according to label per physician discretion, administered as a bolus plus infusion, with the infusion continued for 2-4 hours post-PCI.

SSO₂ Therapy Preparation (Post-PCI)

WARNING: Use aseptic technique throughout procedure.

Only physicians who have received appropriate training in the use of the DownStream System should perform SuperSaturated Oxygen Therapy (SSO₂ Therapy).

SSO₂ Therapy must be performed in the cardiac catheterization laboratory (cath lab). The patient must not be moved while the catheter is in place.

The DownStream System may not be operated concomitantly with magnetic resonance imaging.

The following devices and conditions should be in place prior to starting SSO₂ Therapy:

- Use recommended guidewire to place SSO₂ delivery catheter.
- For a coaxial (single access site) femoral approach (**Figure 4**), place the 7F Merit Medical Custom Sheath Introducer or an 8F introducer sheath with sidearm in the femoral artery.
- For a dual femoral access site, a 5F introducer sheath placed in the contralateral femoral artery is required for blood withdrawal.
- For radial access site approach (**Figure 5**), place SSO₂ delivery catheter through radial sheath introducer. An additional 5F introducer sheath in femoral artery is required for blood withdrawal.

WARNING: Ipsilateral insertion of a second sheath in a single femoral artery for SuperSaturated Oxygen Therapy is strictly contraindicated.

Precaution: The index PCI procedure may be done using femoral or radial access; the SSO₂ Therapy blood draw requires a femoral arterial access site.

Precaution: If the recommended 7F Merit Medical Custom Sheath Introducer is not used, or less than an 8F introducer sheath (coaxial approach) or 5F introducer sheath (dual site approach) is used for femoral access, SSO₂ Therapy may stop due to inadequate flow rate.

- Post PCI draw an arterial blood sample and analyze pO₂. Systemic arterial pO₂ must be greater than or equal to 10.7 kPa or 80 mmHg to proceed with SSO₂ Therapy.

Prepare the Console for Use

1. Position the console and immobilize the wheels by pressing down on and securing the wheel locks. Ensure that the console is plugged into an electrical outlet.
2. Locate the IV pole on the top of the console. Loosen the nut of the IV Pole and slowly raise it to the desired height or until it stops. Tighten the nut to keep it in position. Hang a one-liter bag of sterile physiologic 0.9% normal saline on the console IV pole.
3. To energize the console, rotate the green power lever handle from the Standby position until it latches into the ON position (see **Figure 6**). Engaging the power lever will simultaneously open the oxygen bottle located at the back of the console.

Figure 6. Energizing the DownStream® System Console

4. Confirm that the oxygen bottle contains at least 800 psig (54.4 bar) of oxygen.

Precaution: Do not use the console if the oxygen bottle pressure is less than 800 psig (54.4 bar).

After the console is turned on, the console will start the display. When completed (this may take up to one minute), the display will change to a splash screen (**Figure 7**), indicating the console is ready for use.

Figure 7. Splash Screen

After approximately 15 seconds the screen shown in **Figure 8** appears:


Figure 8. Screen Displays Items Required to Begin Therapy

Load the Cartridge

WARNING: Inspect the DownStream Cartridge and packaging for damage prior to use. Do not use if the cartridge or packaging is damaged.

Prior to use, all packaging, cartridge and lines should be examined for damage. Do not use any defective equipment. Contact TherOx Customer Service for replacement.

1. Open the package. Pass the cartridge tray with the tube set to the sterile field.
2. A user in the sterile field removes the cartridge from the tray, keeps the sterile pouch with tube set and passes cartridge to non-sterile field.
3. Press the button identified as 'CONTINUE' on the display to continue to the next screen and begin the loading process. Once the "CONTINUE" button is pressed the screen shown in **Figure 9** appears:

Figure 9. Screen Displays Prompt to Load Cartridge

4. Open the cartridge compartment door. A small BLUE flashing light in the compartment door indicates the lever for opening the door. When the door is opened, the screen will display the message, "Ready for New Cartridge".

Figure 10. DownStream Cartridge Compartment

5. Align the IV line in the open groove above the cartridge compartment (see **Figure 10**). Ensure the draw and return lines are placed through open grooves in the cartridge compartment.
6. Insert the cartridge, ensuring proper alignment, then close and latch the compartment door. Verify IV tube is not restricted when door is closed. Once the cartridge has been loaded and the door closed, the screen shown in **Figure 11** appears:

Note: The blue light stops flashing.


Figure 11. Prep Screen

7. Insert the IV spike into the bag of sterile physiologic 0.9% normal saline.

WARNING: Use only 0.9% Normal saline (Isotonic solution).

8. Plug the transducer cable into the transducer port.

Precaution: Do not expose electrical connections, including the transducer port, to fluid contact.

9. Press the button identified as 'PREP' on the display to continue to the next screen and begin the cartridge prepping process. Once the 'PREP' button is pressed the screen shown in **Figure 12** appears:

Figure 12. Prep Screen

10. Install the lines through the blood pump, flow probe, return safety clamp, and return line guide while the console is in PREP (see **Figures 13 and 14**).
 - a. Open the blood pump by rotating the blood pump lever upward.
 - b. Insert the draw line into the blood pump, aligning it with the tubing guides (V slots). After achieving proper alignment of the draw line, rotate the blood pump lever downward to close the blood pump.

Note: The blue collar on the draw line must be positioned below the pump head.

- c. Press the release button for the flow probe door. The flow probe door will spring open. Insert the return line in the flow probe slot.

Precaution: Do not stretch the return line while loading into the flow probe.

Precaution: Ensure that the return line is fully seated in the flow probe slot before closing to prevent the flow probe door from pinching the line.

- d. Close the flow probe door.
 - e. Insert the return line through both the return safety clamp and the return line guide.

Note: The red collars should be above and below the return line guide when properly loaded. (**Figure 14**).

Figure 13. Line Installation


Figure 14. Line Installation – Close-up View

11. After the cartridge has been successfully prepped the screen shown in **Figure 15** appears prompting the operator to confirm the correct setup. If user message prompts to spike saline bag, ensure IV spike is fully inserted into saline bag, then press 'CONTINUE' to complete Prep.

Precaution: Do not unplug the transducer once the cartridge has been prepped.


Figure 15. Proper Setup Confirmation Screen

12. After the setup has been confirmed, press the button identified as ‘CONTINUE’ on the display to continue to the next screen and begin the circuit priming process.

Making Patient Connections

WARNING: Do not use any type of catheter for SSO₂ Therapy other than SSO₂ delivery catheter. Inspect the SSO₂ delivery catheter and packaging for damage and kinking prior to use. Do not use if the catheter or packaging is damaged.

SSO₂ delivery catheter preparation

Prior to use, all packaging and devices should be examined for damage. Do not use any defective equipment. Contact TherOx Customer Service for replacement.

1. Transfer the sterile SSO₂ delivery catheter on to the sterile field.
2. Place the SSO₂ delivery catheter into the ostium of the LMCA per standard cardiac catheterization laboratory practice.
3. Remove the guidewire after satisfactory placement of SSO₂ delivery catheter.

Precaution: After placement, inspect the SSO₂ delivery catheter under fluoroscopic visualization to ensure that it is seated correctly in the LMCA and does not restrict blood flow.

4. Maintain a heparinized saline flush or lock through the SSO₂ delivery catheter until connection to the cartridge has been made and the end of the priming sequence.

Precaution: Ensure fluoroscopic contrast agents delivered through the SSO₂ delivery catheter have been flushed prior to priming circuit; these viscous solutions may result in a flow stoppage after connection to the cartridge return line is made.

Cartridge Line Connections

Coaxial Femoral Approach – refer to **Figure 4**

1. Release the ends of the sterile cartridge lines from the pouch on to the sterile field.
2. Attach the draw connector to the sidearm of the arterial sheath, keeping the stopcock on the sidearm closed until ready to prime the circuit with blood.
3. When ready to begin the Prime sequence, open the stopcock on the sidearm of the sheath to allow blood flow into the draw line.

Radial-Femoral Approach – refer to **Figure 5**

1. Release the ends of the sterile cartridge lines from the pouch on to the sterile field.
2. Attach the draw line connector to the sidearm of the femoral sheath, keeping the stopcock on the sidearm closed until ready to prime the circuit with blood.
3. When ready to begin the Prime sequence, open the stopcock on the sidearm of the sheath to allow blood flow into the draw line.

Priming the Blood Path

WARNING: Verify that all connections are secure and tight prior to initiating prime.

Verify that the return line of the cartridge is NOT CONNECTED to the SSO₂ delivery catheter prior to initiating prime.

Precaution: Ensure that draw and return lines are not kinked or restricted prior to starting prime or at any other time during use.

1. Press and continuously hold the green PRIME button located on the console display to blood prime the circuit. As the green PRIME button is held, blood will begin to flow through the circuit and out of the return line.
2. As the console pump begins to operate, a message appears on the screen indicating that the circuit is priming (see **Figure 16**). The blood pump will slow down to facilitate priming the return line.


Figure 16. Prompt to Press and Hold PRIME Button

3. When the return line is blood primed and bubble free, the physician makes a wet-to-wet connection with the red return line connector to the SSO₂ delivery catheter hub.

WARNING: The bubble detector is disabled during prime. Ensure that return line is fully blood primed with no presence of air or bubbles in the line prior to making wet-to-wet connection to the SSO₂ delivery catheter.

4. Keep pressing the green PRIME button until the screen gives notice that the console pressure is stabilizing (see **Figure 17**) and allows the release of the green PRIME button.


Figure 17. Prompt to Release PRIME Button

5. At the completion of prime, visually inspect that all line connections are leak free.

WARNING: Do not disconnect return line from SSO₂ delivery catheter or draw line from sheath when therapy is in progress. Procedure must be stopped before disconnection.

Start and Monitor SSO₂ Therapy

Note: SSO₂ Therapy is performed in the cardiac catheterization laboratory.

1. At completion of Prime, the following screen will appear (**Figure 18**):

Figure 18. Screen Display for Entering SSO₂ ON

2. Press the button on the display identified as “SSO₂ ON” to begin SSO₂ Therapy.
3. SSO₂ Therapy will continue for 60 minutes. The screen shown in **Figure 19** will be displayed during therapy.

Figure 19. Screen Display During SSO₂ Therapy

Terminate SSO₂ Therapy and Blood Flow

After 60 minutes of treatment, SSO₂ Therapy will stop and the console will indicate (on the display) that the infusion is complete. The blood pump will continue to circulate normoxic arterial blood through the cartridge; the blood will not be mixed with SSO₂ solution (see **Figure 20**).

Figure 20. Screen Display after SSO₂ Therapy Complete

1. Press the button identified as ‘END PROCEDURE’ on the display.
2. The yellow message area on the screen will display a prompt to confirm the user request to end the procedure.
3. Press the button ‘END PROCEDURE’ again to indicate procedure completion and to stop flow through the cartridge.

- The display will change to a screen that prompts the operator to disconnect from the patient and to press the 'CONTINUE' button when finished. See **Figure 21**.

Note: Decompression in the cartridge chamber will be accompanied by an audible hissing sound. Do not attempt to open the cartridge compartment door until prompted.


Figure 21. Screen Display Continue Button

- Once the blood pump has stopped, close the stopcock on the sheath sidearm.
- Detach the SSO₂ delivery catheter from the cartridge return line by disconnecting the red return connector from the SSO₂ delivery catheter hub.
- Disconnect the draw connector from the sheath sidearm.
- Immediately remove the SSO₂ delivery catheter from the patient.

Note: It is not recommended to return the blood in the cartridge to the patient. It is recommended to leave the blood (approximately 60 ml) in the cartridge for disposal (using standard hospital procedures to eliminate hazardous waste).

Unload the Cartridge

- Press the button on the display identified as 'CONTINUE' to unload the cartridge. The display will change to a screen that notifies the operator that the console is preparing to unload the cartridge (see **Figure 22**).


Figure 22. Preparing Cartridge for Unload

- Slide the blood pump lever up toward the top of the console and remove the draw line.
- Open the flow probe door by pressing the flow probe door release button and remove the return line.
- Remove the return line first from the return safety clamp and then from the return line guide.
- When the screen displays the "Open Door and Unload Cartridge" prompt (Figure 23), open the cartridge compartment door, unplug the transducer cable from the transducer port, and then remove the cartridge.


Figure 23. Prompt to Unload Cartridge

6. After the cartridge has been removed, close and secure the compartment door.
7. Dispose of the single-use cartridge and sterile physiologic 0.9% normal saline bag per standard hospital procedure.

WARNING: **Do not re-use the single-use cartridge and SSO₂ delivery catheter.**
 Dispose of used cartridges and SSO₂ delivery catheters using standard hospital procedures to eliminate hazardous waste.

Shut Down the Console

Move the green power lever handle from the ON position to the Standby position. This will also close the oxygen bottle located at the back of the console.

Emergency Stop

Press the large red Emergency Stop switch on the front panel, at any time, to immediately stop console operation (See **Figure 24**). To resume console use, release the Emergency Stop switch by turning it clockwise.

Note: The cartridge must be replaced to resume treatment after an Emergency Stop.

WARNING: **Under no circumstances should the extracorporeal circuit be disconnected and restarted outside the cath lab.**

Figure 24. Emergency Stop Switch

Responding To Error Messages

Under certain conditions, an error message may appear in the message area in the center of the display screen. If the error message appears in yellow, follow the recovery procedures indicated on the screen. If the error message appears in red, the console is stopped. Follow the instructions displayed in the message area on the display.

Help Screen

The Help screen provides user information specific to the current step in the operating sequence. The Help screen also displays the TherOx Technical Support telephone number.

Equipment Maintenance - DownStream System Console

WARNING: Only qualified personnel trained in oxygen safety and handling should change the oxygen bottle.

- As needed, wipe the console cabinet and the interior of the cartridge compartment with an approved hospital bactericidal agent. Avoid fluid contact with the transducer port.
- To maintain batteries in fully charged condition, the console must be powered from an electrical outlet when not in use.
- After a low battery warning, allow at least six hours to fully recharge the batteries before the next use.
- Expected battery life is 200 cycles or two years.
- Consoles that have been in storage (un-powered) should be charged for at least 12 hours before use.
- In normal use, a full oxygen bottle should provide over 50 uses. When the oxygen bottle contains less than 800 psig (54.4 bar), it should be replaced by qualified personnel per standard hospital procedures.
- The console contains no hospital serviceable parts.

Installation and Service Requirements

Only TherOx personnel or representatives are authorized to install and service the DownStream Console.

To ensure safe operation of the TherOx console, preventive maintenance including a comprehensive technical inspection and electrical safety test is required annually.

Device Labeling Symbols

	Contains or presence of: Bis(2-ethylhexyl) phthalate (DEHP)
	Catalog Number
	Read all instructions prior to use
	Batch Code
	Sterilized with ethylene oxide gas
	Use By (YYYY-MM-DD)
	Do Not Re-use
	Storage Temperature
	Keep dry
	Cardiovascular floating defibrillation proof
	A/C Current
	Amps
	Volts
	Power is on
	Console is in standby mode, and the battery is charging (while plugged into an electrical outlet)
	Fusing
	Electrostatic sensitive devices

Refer to instruction manual/booklet

Equipotentiality

Manufacturer


Oxygen Connector

Maximum Oxygen inlet pressure

Maximum Console weight

Date of manufacture (YYYY-MM-DD)

Do not use if package is damaged

Consult the instructions for use for important cautionary information

Separate collection for waste electrical and electronic equipment (WEEE)

If you wish to discard electrical and electronic equipment (EEE), please contact your distributor for further information.

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TherOx, Inc.